

TG TUBE FURNACE

GENERAL INFORMATION

The furnace operates at temperatures up to 2500°C in an inert or vacuum atmosphere. The insulation and shielding is formed with pressed graphite compounds to withstand high temperatures.

The furnace is typically custom-built as per customer request. Please request a quote through the Mellen Sales Team for a specific furnace.

HEATING ELEMENT

The furnace utilizes Graphite heating elements capable of 2500°C in an inert atmosphere. These heating elements are for high temperature applications and can be custom-fit to a specific size and orientation.

CONTROLLER

Sold Separately.



SERIES	TG Tube Furnace
MAX TEMPERATURE	2500°C
ELEMENT TYPE	Graphite
CONTROLLER	Sold Separately
VOLTS	120V & 208V/240V

*As always, if this product does not meet your requirements, please contact us.
Mellen will work with you personally to create a product that best suits your needs every time.*

Prices are valid until 2018.

