

TI TUBE FURNACE

GENERAL INFORMATION


The furnace operates at temperatures up to 1450°C. The furnace pictured has an ultra-small 1 inch zone, making it ideal for miniature applications, fine areas, extreme gradients, lab services, and zone refining. The TI Tube Furnace is ideal for a laboratory environment where footprint is critical.

HEATING ELEMENT

The furnace utilizes Sunfire heating elements capable of up to 1450°C. These heating elements are coiled and spread to best radiate heat in a confined area.

CONTROLLER

Sold Separately.



SERIES	TI Tube Furnace
MAX TEMPERATURE	1450°C
ELEMENT TYPE	Sunfire
CONTROLLER	Sold Separately
VOLTS	120V & 208V/240V

*As always, if this product does not meet your requirements, please contact us.
Mellen will work with you personally to create a product that best suits your needs every time.*

Prices are valid until 2018.

